

bionet

Veterinary wearable device that can wirelessly measure pet vital signs

O'Pet Curv



Veterinary Wearable Smart Patch O'Pet Curv can wirelessly monitor pet vital signs such as heart rate, respiratory rate, and body temperature through an app that comes with it.

When the measured vital signs are outside the set range, an alarm goes off to act immediately, and the measured values are stored on the cloud server in real time to check the patient's condition regardless of time and place.

Now save time and money to monitor high-risk patients with O'Pet Curv.

Never miss the critical moment to save a pet's life.



Wireless
Wearable
Sensor

LIVE

Heart Rate
Respiratory Rate
Skin Temperature



Alarm



Remote
Patient
Monitoring



Vital Sign
Chart
ECG Record



Battery
72 Hours



IP68
Waterproof

Specification



O'Pet Curv
(Wireless Wearable
Sensor)



OA-100
(Hydrogel
Adhesive)




OC-100
(Charger)



OV-101
Reusable electrode



OV-102
Holding strap




Size, Weight	90 x 34 x 8.5 mm, 18g
ECG	1 Channel, 125 Hz
Heart Rate	40 ~ 300 BPM
Respiration Rate	5 ~ 60 RPM
Skin Temperature	10 ~ 45 °C
Ingress Level	IP68
Battery	4 hour charging, 72 hour operation
Communication	Bluetooth 5.1
Software	App for Android/iOS/Web Browser

Overall Usage



Pairing O'Pet Curv
to mobile device




Attach a hydrogel
adhesive to O'Pet Curv



Attach O'Pet Curv
to the chest
of companion animal



Display and record
the vital signs.




Review the recorded
vital signs and ECG

Monitoring Display, Vital Sign Chart, ECG Chart



Remote Monitoring Platform



Inside Hospital



Recovery Room
or Operation Room

Remote Central
Monitoring Room



O'Pet LinkVue



- Detects heart rate, respiratory rate, and skin temperature from biosignals acquired through the O'Pet Curv and displays them on the screen.
- An alarm automatically goes off if the reading is outside the set alarm range.
- Saves vital signs and ECG to the mobile device.
- Connects and monitors up to 6 wireless sensors at the same time.
- Monitors the patient remotely by “monitoring sharing” function.

O'Pet StationVue



- Remotely monitors vital signs of all wireless devices connected to O'Pet LinkVue.
- Monitors patients' conditions and receives alarms in real time through PCs and mobile devices inside and outside the hospital.
- When used on a PC, the patient's vital sign and ECG records can be downloaded and stored on your PC to attach to an EMR chart.
- Operates on Web browsers with <https://sv.openbio.io>


Applications



**Emergency patient monitoring
with no movement**



**Monitoring patients under anesthesia
during surgery**



**Monitoring critical patients recovering
from surgery**

Advantage & Benefit

01. Wireless wearable sensors make it quick and easy to apply to patients.

- Even in emergency situations, sensors can be quickly and easily attached to the patient's body, and the patient does not reject or feel uncomfortable with the attached sensor.
- You can reliably monitor vital signs from aggressive or irritating patients without being exposed to dangerous situations.
- There is no wire and the sensor is attached to the patient's body, providing excellent quality vital signs.

02. The alarm function allows you to respond quickly.

- An alarm is triggered on an app running on your mobile device or PC when the reading is outside the set alarm range.
- The alarm function can be used to respond quickly before the patient's condition gets worse.
- When an alarm is triggered above the heart rate, you can set up a 40-second ECG to be automatically saved on the cloud server.

03. Remote monitoring function allows you to check the patient's vital signs anytime, anywhere.

- Cloud-based apps allow you to monitor your patient's condition anytime, anywhere, and review and share the saved vital sign trends and ECG records.

04. Lower hospital operating costs and increase work efficiency.

- O'Pet monitors motionless patients continuously, so you can focus on other tasks with confidence.
- By automatically recording vital signs, you can save money on measuring and recording vital signs.

bionet



Headquarters Bionet Co., Ltd.

5F, Dream Markone 61 Digital-ro 3I-gil, Guro-gu, Seoul, Korea 08375

TEL. +82 2 6292 6410

E-Mail. info.sales@ebionet.com

www.ebionet.com